



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Mouse LE/AF Purified Anti-Mouse IL-12/IL-23 p40 Monoclonal Antibody

单克隆抗体

产品货号: MA02423

第三版

免疫原:None

亚型:Rat IgG2a

克隆:C17.8

描述:The C17.8 antibody reacts with the free p40 monomer of IL-12 as well as the p40 subunit complexed as homodimer or complexed with p35 as a p70 heterodimer. Clone C17.8 has also been shown to be cross-reactive with mouse IL-23, which also contains the p40 subunit. The immunogen used to generate C17.8 hybridoma was recombinant mouse IL-12 p70.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF