

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Rat LE/AF Purified Anti-Rat CD4 (domain 1) Monoclonal Antibody

单克隆抗体

产品货号: MA02427

第三版

免疫原:None

亚型:Mouse IgG2a

克隆:OX-38

描述:CD4, also known as T4, is a 55kD glycoprotein member of the immunoglobulin superfamily and is expressed on majority of thymocytes, macrophages, and a peripheral T cell subset (T helper cells). CD4 is a T cell co-receptor that interacts with the MHC class II molecule and is involved in T cell activation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF