

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat AGTR1 Monoclonal Antibody

单克隆抗体

产品货号: MA02429

第三版

免疫原:Fusion Protein of AGTR1

亚型:IgG2a

克隆:2C5G11

**描述:**Angiotensin II (Ang II),the main effector molecule of the renin-angiotensin system,exerts its actions mainly via interaction with type-1 angiotensin II receptor (AGTR1,also named as AT1R),thereby contributing to blood pressure regulation. AGTR1 mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system. By regulating vascular tone,cardiovascular function,salt and water homeostasis,AGTR1 exerts an indispensable physiological role . AGTR1 has been implicated in diverse aspects of human disease,from the regulation of blood pressure and cardiovascular homeostasis to cancer progression .

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF