



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Pig ALOX5 Monoclonal Antibody

单克隆抗体

产品货号: MA02430

第三版

免疫原:Fusion Protein of ALOX5

亚型:IgG2a

克隆:3B1B7

描述:ALOX5(arachidonate 5-lipoxygenase) is also named as LOG5,5-LO and belongs to the lipoxygenase family. It catalyzes the first step in leukotriene biosynthesis,and plays a role in inflammatory processes. It is a 78 kDa iron-containing enzyme that translocates from the cytosol to associate with FLAP,and converts arachidonate in two steps into the unstable intermediate LTA4. It is essential for cysteinyl-leukotriene (cys-LT) production,critical mediators in asthma.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF