

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse ANXA1 Monoclonal Antibody

单克隆抗体

产品货号: MA02432

第三版

免疫原:Fusion protein of Annexin A1

亚型:IgG1

克隆:1E1B7

描述:ANXA1,also named as ANX1,LPC1 and p35,belongs to the annexin family. It is a member of the annexin family of Ca²⁺-binding and phospholipidbinding proteins,and it is particularly abundant in various populations of peripheral blood leukocytes. ANXA1 regulates phospholipase A2 activity. It seems to bind from two to four calcium ions with high affinity. ANXA1 is a pleiotrophic protein produced by many cell types including peripheral blood leucocytes.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF