



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ANXA3 Monoclonal Antibody

单克隆抗体

产品货号: MA02433

第三版

免疫原:Fusion protein of Annexin A3

亚型:IgG2a

克隆:1A11E4

描述:Annexin A3 (Anxa3),also named as lipocortin3,placental anticoagulant protein3 (PAP-III),belongs to annexin family. Annexins are calcium-dependent phospholipids-binding high-abundant intracellular proteins. ANXA3 has been reported to play crucial roles in a diverse range of molecular and cellular processes,including ion channel regulation,anticoagulation,anti-inflammation,vesicle traffick-ing,signal transduction,cell migration,and cell proliferation. It also exhibits important roles in cancer progression,metastasis,and drug resistance.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF