

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BCL6 Monoclonal Antibody

单克隆抗体

产品货号: MA02437

第三版

免疫原:Fusion protein of BCL6

亚型:IgG1

克隆:1E6B1

描述:BCL6,a zinc finger transcription factor,contains a N-terminal BTB/POZ domain and C-terminal zinc finger DNA-binding motifs and represses transcription of a wide range of target proteins and microRNAs. BCL6 protein has been reported as a master regulator of B lymphocyte development and growth,and altered BCL6 protein expression was implicated in pathogenesis of diverse human hematologic malignancies,especially in the diffuse large B cell lymphoma (DLBCL). BCL6 is required for the development of T follicular helper T cells (TFH),a helper T cell subset required for the formation of mature and productive GCs. BCL6 has also been shown to play important regulatory roles in macrophages.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF