

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human BMP5 Monoclonal Antibody

单克隆抗体

产品货号: MA02439

第三版

免疫原:Fusion Protein of BMP5

亚型:IgG1

克隆:1G6C10

描述:Bone morphogenetic protein-5 (BMP5),is a member of the transforming growth factor-beta (TGF-beta) superfamily. The superfamily includes large families of growth and differentiation factors. Bone morphogenetic proteins were originally identified by an ability of demineralized bone extract to induce endochondral osteogenesis in vivo in an extraskelatal site. These proteins are synthesized as prepropeptides,cleaved,and then processed into dimeric proteins. BMP5 may act as an important signaling molecule within the trabecular meshwork and optic nerve head,and may play a potential role in glaucoma pathogenesis. BMP5 is differentially regulated in multiple human cancers.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF