



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat BABAM2 Monoclonal Antibody

单克隆抗体

产品货号: MA02440

第三版

免疫原:Fusion Protein of BRE

亚型:IgG2a

克隆:2E5A5

描述:BRE,also named as BRCC45,is a death receptor-associated antiapoptotic protein,which is also involved in DNA-damage repair,and K63-specific deubiquitination. BRE and its XIAP-sustaining mechanism could represent novel targets for anti-cancer therapy. BRE may regulate TNF-alpha signaling through its interactions with TNFRSF1A. It plays a role in homeostasis or cellular differentiation in cells of neural,epithelial and germline origins. BRE gene has four isoforms with MW 42-47 kDa.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF