



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CALB1 Monoclonal Antibody

单克隆抗体

产品货号: MA02441

第三版

免疫原:Fusion Protein of Calbindin-D28k

亚型:IgG1

克隆:1F8B9

描述:Calbindin-D28k is a member of the calcium-binding protein superfamily that includes calmodulin and troponin C. It is a cytosolic calcium binding protein highly expressed in the distal tubule,intestines,central nervous system,primary murine osteoblast cells and in several other organs. It plays an important role in the intracellular calcium homeostasis,its strong buffering capacity prevents cytotoxic effect of high concentration of free calcium in kidney,brain,pancreas,and intestine.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF