



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CALCOCO2 Monoclonal Antibody

单克隆抗体

产品货号: MA02442

第三版

免疫原:Fusion Protein of CALCOCO2

亚型:IgG2a

克隆:2F2C2

描述:CALCOCO2,also named as NDP52,is an autophagy receptor. It plays a role in ruffle formation and actin cytoskeleton organization. Mouse/Rat Calcoco2 has some isoforms with MW 28-40 kDa and 67 kDa. Human CALCOCO2 has some isoforms with MW 43-47 kDa and 52-55 kDa,

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF