



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CAMSAP2 Monoclonal Antibody

单克隆抗体

产品货号: MA02444

第三版

免疫原:Fusion Protein of CAMSAP2

亚型:IgG1

克隆:2A10A1

描述:CAMSAP2,also known as CAMSAP1L1,is a microtubule minus-end binding protein,and is required for proper microtubule organization in neurons. CAMSAP2 stabilizes noncentrosomal microtubules and is important for neuronal polarity,axon specification,and dendritic branch formation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF