



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD2 Monoclonal Antibody

单克隆抗体

产品货号: MA02451

第三版

免疫原:Fusion Protein of CD2

亚型:IgG1

克隆:3A10B2

描述:CD2 is a cell surface glycoprotein present on a majority of thymocytes,all mature T cells and subset of NK cells but not on B lymphocytes. It is a pan T-cell marker. CD2 interacts with lymphocyte function-associated antigen (LFA-3/CD58) and CD48/BCM1 to mediate adhesion between T-cells and other cell types. CD2 is implicated in the triggering of T-cells,the cytoplasmic domain is implicated in the signaling function.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF