

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CCP110 Monoclonal Antibody

单克隆抗体

产品货号: MA02460

第三版

免疫原:Fusion Protein of CP110

亚型:IgG2a

克隆:1A3C12

**描述:**CP110,also named as CCP110 and KIAA0419,is an 110kDa protein. This gene is different to the gene CEP110 (geneID:11064; CNTRL). CP110 is a centriolar protein that positively regulates centriole duplication while restricting centriole elongation and ciliogenesis. And it acts as a key negative regulator of ciliogenesis in collaboration with CEP97 by capping the mother centriole thereby preventing cilia formation.

**配方:**

**浓度:**None

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF