



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Cyclin A2 Monoclonal Antibody

单克隆抗体

产品货号: MA02463

第三版

免疫原:Fusion Protein of Cyclin A2

亚型:IgG2a

克隆:1C6C11

描述:CCNA2 belongs to the highly conserved cyclin family, whose members are characterized by a dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases. Different cyclins exhibit distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event. In contrast to cyclin A1, which is present only in germ cells, this cyclin is expressed in all tissues tested. CCNA2 binds and activates CDC2 or CDK2 kinases, and thus promotes both cell cycle G1/S and G2/M transitions.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF