

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat DPT Monoclonal Antibody

单克隆抗体

产品货号: MA02467

第三版

免疫原:Fusion Protein of Dermato pontin

亚型:IgG2b

克隆:4B1A12

描述:DPT (dermatopontin) is an extracellular matrix protein belonging to the dermatopontin family. Skin appears to be the richest source of dermatopontin, comprising about 15 mg/kg of the wet weight. Expression of dermatopontin is not limited to connective tissue, as investigations show that dermatopontin mRNA is expressed in cultured fibroblasts, muscle, heart, pancreas, and lung. Dermato pontin comprises a considerable proportion of the noncollagenous extracellular matrix proteins. It is involved in promotion of cellular adhesion, extracellular matrix assembly, wound healing, and positive modification of the growth inhibition activity of TGF- β 1.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF