

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human EPCAM Monoclonal Antibody

单克隆抗体

产品货号: MA02471

第三版

免疫原:Fusion Protein of EPCAM

亚型:IgG2a

克隆:2A2D5

**描述:**Epithelial cell adhesion molecule (EpCAM,CD326) is a type I transmembrane glycoprotein that functions as a homophilic,epithelial-specific intercellular cell-adhesion molecule. In addition to cell adhesion,EpCAM is also involved in cellular signaling,cell migration,proliferation,and differentiation. EpCAM is highly expressed on most carcinomas and therefore of potential use as a diagnostic and prognostic marker for a variety of carcinomas,and has become a therapeutic target.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF