



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse FOXO3 Monoclonal Antibody

单克隆抗体

产品货号: MA02476

第三版

免疫原:Fusion Protein of FOXO3A

亚型:IgG2a

克隆:1F12D11

**描述:**In mammalian cells,FoxO factors act as potent transcriptional activators that upregulate the expression of programs of genes involved in stress resistance,cell cycle arrest,differentiation,apoptosis,autophagy,and metabolism. FOXO3 has a key role in the maintenance of adult neural and hematopoietic stem cells. It also associates with extreme longevity in humans and regulates the homeostasis of adult stem cell pools in mammals,which may contribute to longevity. Phosphorylation inhibits FOXO3 nuclear translocation that is essential for its regulation and function

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF