

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat GDF10 Monoclonal Antibody

单克隆抗体

产品货号: MA02480

第三版

免疫原:Fusion Protein of GDF10

亚型:IgG2a

克隆:1B7D12

描述:Growth differentiation factor-10 (GDF10),also known as BMP3b,is a member of the transforming growth factor-b (TGF-b) superfamily. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer. This promotes neural repair after stroke. GDF10 plays a role in head formation and may have multiple roles in skeletal morphogenesis. Additionally,GDF10 may act as a tumor suppressor and reduced expression of this gene is associated with oral cancer.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF