



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GNPAT Monoclonal Antibody

单克隆抗体

产品货号: MA02483

第三版

免疫原:Fusion Protein of GNPAT

亚型:IgG2b

克隆:2E5C10

描述:GNPAT,also named as DAPAT and DHAPAT,belongs to the GPAT/DAPAT family. It is a key enzyme in the biosynthesis of ether phospholipids. GNPAT is localized exclusively within peroxisomes. Full GNPAT activity depends not only on the presence of AGPS,but also on the integrity of substrate channeling from GNPAT to AGPS. This antibody recognize the 65-69kd GNPAT protein.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF