

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human GOLM1 Monoclonal Antibody

单克隆抗体

产品货号: MA02484

第三版

免疫原:Fusion Protein of GP73/GOLPH2

亚型:IgG1

克隆:3D1F1

描述:Golgi phosphoprotein 2 (GOLPH2,also known as GP73 or GOLM1) is a resident Golgi type-II membrane protein. It is predominantly expressed in the epithelial cells of many human tissues. GOLPH2 traffics through endosomes and can be secreted into the circulation. Its expression is upregulated in a number of tumors and GOLPH2 could be a promising serum marker for hepatocellular carcinoma.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF