

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse CD101 Monoclonal Antibody

单克隆抗体

产品货号: MA02489

第三版

免疫原:Fusion Protein of IGSF2

亚型:IgG1

克隆:1E6A7

描述:IGSF2,also known as CD101,is a low abundance protein. Expression of CD101 is found on monocytes,granulocytes and dendritic cells. In addition expression on T lymphocytes is important for cell activation. It plays a role as inhibitor of T-cells proliferation induced by CD3 and inhibits expression of IL2RA on activated T-cells and secretion of IL2.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF