



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL17A Monoclonal Antibody

单克隆抗体

产品货号: MA02492

第三版

免疫原:Fusion Protein of IL17a

亚型:IgG2b

克隆:2F9D1

描述:IL17a is a pro-inflammatory cytokine as a member of the interleukin-17 family. It plays a central role in host defense against diverse pathogens. It is produced by activated T-cells and certain cell types of innate immune system. The active protein functions as either a homodimer with other interleukin-17 family members and signals through the interleukin-17 receptor to induce inflammatory cytokine production. Aberrant expression of this gene is associated with autoinflammatory diseases including rheumatoid arthritis, psoriasis and multiple sclerosis.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF