

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat IRF4 Monoclonal Antibody

单克隆抗体

产品货号: MA02494

第三版

免疫原:Fusion Protein of IRF4

亚型:IgG1

克隆:1A12D5

描述:Interferon regulatory factor 4 (IRF4) is a member of the IRF family of transcription factors,expressed in most cell types of the immune system. IRF4 is a transcription factor essential for the development of T helper-2 (Th2) cells,IL17-producing Th17 cells,and IL9-producing Th9 cells,as well as dendritic cell (DC). It binds to and activates the interferon-stimulated response element (ISRE) of the MHC class I promoter,and may have a role in ISRE-targeted signal transduction mechanisms specific to lymphoid cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF