



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat LDHB Monoclonal Antibody

单克隆抗体

产品货号: MA02499

第三版

免疫原:Fusion Protein of LDHB

亚型:IgG1

克隆:5E9G1

描述:LDHB,also named as LDH-H and NY-REN-46,belongs to the LDH/MDH superfamily. And LDH family. It is an enzyme which catalyzes the reversible conversion of lactate and pyruvate,and NAD and NADH,in the glycolytic pathway.It is glycolytic enzymes and it is positively regulated by mTOR.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF