



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat MEA1 Monoclonal Antibody

单克隆抗体

产品货号: MA02503

第三版

免疫原:Fusion Protein of MEA1

亚型:IgG1

克隆:

描述:MEA1,also named as Male-enhanced antigen 1,is a 185 amino acid protein,which is highly expressed in testis. MEA1 may play an important role in spermatogenesis and/or testis development.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF