



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MKP-2 Monoclonal Antibody

单克隆抗体

产品货号: MA02506

第三版

免疫原:Fusion Protein of MKP-2

亚型:IgG1

克隆:1G1C12

描述:DUSP4,also named as MKP-2,is a member of the mitogen-activated protein (MAP) kinase phosphatase family which has been suggested to play an important role in the feedback control of MAP kinase-mediated gene expression. DUSP4 preferentially inactivates ERK1,ERK2 and JNK,is expressed in a variety of tissues,and is localized in the nucleus.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF