

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MSLN Monoclonal Antibody

单克隆抗体

产品货号: MA02509

第三版

免疫原:Fusion Protein of MSLN

亚型:IgG1

克隆:2E3A4

**描述:**This gene encodes a precursor protein that is cleaved into two products, megakaryocyte potentiating factor and mesothelin. Megakaryocyte potentiation factor functions as a cytokine that can stimulate colony formation in bone marrow megakaryocytes. Mesothelin is a glycosylphosphatidylinositol-anchored cell-surface protein that may function as a cell adhesion protein. This protein is overexpressed in epithelial mesotheliomas, ovarian cancers and in specific squamous cell carcinomas. Alternative splicing results in multiple transcript variants.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF