

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MUM1 Monoclonal Antibody

单克隆抗体

产品货号: MA02510

第三版

免疫原:Fusion Protein of MUM1

亚型:IgG2a

克隆:3B1D2

描述:MUM1,also named as EXPAND1,is an architectural component of the chromatin,which in response to DNA damage serves as an accessory factor to promote cell survival. It is a 53BP1-associated DNA damage-responsive factor. MUM1/EXPAND1 is different to another protein IRF4/MUM1.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF