



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MYOC Monoclonal Antibody

单克隆抗体

产品货号: MA02511

第三版

免疫原:Fusion Protein of MYOC

亚型:IgG2a

克隆:4B1F6

描述:Myocilin (MYOC),also known as TIGR or GLC1A,is a secreted glycoprotein of the olfactomedin family,originally identified in trabecular meshwork cells after prolonged treatment with glucocorticoids,and,independently,in the retina . It may participate in the obstruction of fluid outflow in the trabecular meshwork. Defects in MYOC are the cause of primary open angle glaucoma (POAG) .

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF