

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat NONO Monoclonal Antibody

单克隆抗体

产品货号: MA02516

第三版

免疫原:Fusion Protein of NONO

亚型:IgG1

克隆:2A2B10

描述:Non-pou domain-containing octamer-binding protein(NONO),other named p54nrb,belongs to DBHS family,whose protein are predominately nuclear and involves in transcriptional and posttranscriptional gene regulatory functions as well as DNA repair. This family has roles in a wide range of biological processes,such as the regulation of circadian rhythm,carcinogenesis and progression of cancer. NONO forms heterodimer with SFPQ,and has been implicated in nuclear retention of A- to I-edited RNA as paraspeckle components,pre-mRAN 3'-end formation,cAMP cycling and transcriptional activation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF