
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse NUDT21 Monoclonal Antibody

单克隆抗体

产品货号: MA02519

第三版

免疫原:Fusion Protein of NUDT21

亚型:IgG1

克隆:1C4C8

描述:Cleavage and polyadenylation specific factor 5,25 kDa (CPSF5,synonym: CFIM25) is one subunit of a cleavage factor required for 3' RNA cleavage and polyadenylation processing. The interaction of CPSF5 with the RNA is one of the earliest steps in the assembly of the 3' end processing complex and facilitates the recruitment of other processing factors.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF