



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat NUP155 Monoclonal Antibody

单克隆抗体

产品货号: MA02520

第三版

免疫原:Fusion Protein of NUP155

亚型:IgG2a

克隆:1A7F2

描述:NUP155 is a nucleoporin,the major component of the nuclear pore complex (NPC) that is essential for nuclear envelope biogenesis. Mutation of NUP155 has been linked to cardiovascular disease.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF