

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human OGN Monoclonal Antibody

单克隆抗体

产品货号: MA02522

第三版

免疫原:Fusion Protein of OGN

亚型:IgG1

克隆:2A5F11

描述:This gene encodes a member of the small leucine-rich proteoglycan (SLRP) family of proteins. The encoded protein induces ectopic bone formation in conjunction with transforming growth factor beta and may regulate osteoblast differentiation. High expression of the encoded protein may be associated with elevated heart left ventricular mass. Alternative splicing results in multiple transcript variants.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF