

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human OIP5 Monoclonal Antibody

单克隆抗体

产品货号: MA02523

第三版

免疫原:Fusion Protein of OIP5

亚型:IgG2b

克隆:1E9B5

描述:The protein encoded by this gene localizes to centromeres, where it is essential for recruitment of CENP-A through the mediator Holliday junction recognition protein. Expression of this gene is upregulated in several cancers, making it a putative therapeutic target. Two transcript variants encoding different isoforms have been found for this gene.OIP5 (Opa Interacting Protein 5) is a Protein Coding gene. Among its related pathways are Chromosome Maintenance and Cell Cycle, Mitotic.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF