



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat P2RX4 Monoclonal Antibody

单克隆抗体

产品货号: MA02524

第三版

免疫原:Fusion Protein of P2RX4

亚型:IgG2b

克隆:1A5A6

描述:Nucleotides,the structural subunits of the nucleic acids,are also important extracellular signaling molecules. P2 receptors are a family of cell surface receptors that mediate a wide variety of physiologic effects in response to extracellular nucleotides. These receptors fall into two classes: P2X receptors,which are ligand-gated ion channels that mediate calcium and potassium fluxes in response to ATP,and P2Y receptors,which are G protein-coupled receptors (GPCRs). P2RX4 (P2X purinoceptor 4) is a receptor for ATP acting as a ligand gated ion channel.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF