

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat P2RX5 Monoclonal Antibody

单克隆抗体

产品货号: MA02525

第三版

免疫原:Fusion Protein of P2RX5

亚型:IgG2b

克隆:2A2A6

**描述:**The product of this gene belongs to the family of purinoceptors for ATP. This receptor functions as a ligand-gated ion channel. Alternative splicing results in multiple transcript variants. Read-through transcription also exists between this gene and the neighboring downstream gene, TAX1BP3.P2RX5 (Purinergic Receptor P2X 5) is a Protein Coding gene. Diseases associated with P2RX5 include Cystinosis. Among its related pathways are Calcium signaling pathway and CREB Pathway. GO annotations related to this gene include transmembrane signaling receptor activity and purinergic nucleotide receptor activity. An important paralog of this gene is P2RX5-TAX1BP3.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF