



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PADI2 Monoclonal Antibody

单克隆抗体

产品货号: MA02526

第三版

免疫原:Fusion Protein of PADI2

亚型:IgG2b

克隆:2D8C12

描述:PADI2,also named as KIAA0994,PDI2,PAD-H19 and PAD2(Peptidylarginine deiminase II ),belongs to the protein arginine deiminase family. It catalyzes the deimination of arginine residues of proteins. PADI2 may play a regulatory role in the expression of lactation related genes via histone citrullination during diestrus . PADI2 has two isoforms with MW 75 kDa and 49 kDa.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF