

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PDCD10 Monoclonal Antibody

单克隆抗体

产品货号: MA02527

第三版

免疫原:Fusion Protein of PDCD10

亚型:IgG1

克隆:1E5A8

描述:PDCD10,also named CCM3 and TFAR15,belongs to the PDCD10 family. PDCD10 promotes cell proliferation,increases MAPK activity and modulates apoptotic pathways. PDCD10 is important for cell migration and for normal structure and assembly of the Golgi complex. PDCD10 is required for normal angiogenesis,vasculogenesis and hematopoiesis during embryonic development,moreover,defects in PDCD10 showed abnormal cardiovascular development.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF