



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PEBP1 Monoclonal Antibody

单克隆抗体

产品货号: MA02530

第三版

免疫原:Fusion Protein of PEBP1

亚型:IgG1

克隆:2D4D11

描述:Raf kinase inhibitor protein (RKIP),also termed phosphatidylethanolamine binding protein PEBP1,was initially identified as a potent inhibitor of Raf-1/MEK/ERK,NF- κ B,and G-protein-coupled receptor signaling pathways. Later RKIP has been further identified as a metastasis suppressor and its loss of expression has been reported in various cancers. Its expression has been proposed as a prognostic marker for patients diagnosed with the above cancers.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF