



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat POSTN Monoclonal Antibody

单克隆抗体

产品货号: MA02531

第三版

免疫原:Fusion Protein of Periostin

亚型:IgG1

克隆:1A11A3

描述:Periostin (POSTN,PN),originally named as osteoblast-specific factor 2 (OSF-2),is a 90-kDa secreted protein which is now classified as a matricellular protein. It is present in a wide variety of normal adult tissues and fetal tissues,and has a role in bone,tooth and heart development and function. Studies show that periostin is overexpressed in a broad range of human cancer types,includidng lung,ovary,breast and colon cancers. Recent evidence reveals that periostin is expressed by fibroblasts in the normal tissue and in the stroma of the primary tumour,and it is required to allow cancer stem cell maintenance.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF