

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PPARGC1A Monoclonal Antibody

单克隆抗体

产品货号: MA02542

第三版

免疫原:Fusion Protein of PPARGC1A

亚型:IgG1

克隆:1C1B2

描述:PPARGC1A,also named as Peroxisome proliferator-activated receptor gamma coactivator 1-alpha,is a 798 amino acid protein,which Contains 1 RRM (RNA recognition motif) domain and localizes in the nucleus. PPARGC1A is transcriptional coactivator for steroid receptors and nuclear receptors. PPARGC1A can regulate key mitochondrial genes that contribute to the program of adaptive thermogenesis and plays an essential role in metabolic reprogramming in response to dietary availability through coordination of the expression of a wide array of genes involved in glucose and fatty acid metabolism.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF