

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RAB5A Monoclonal Antibody

单克隆抗体

产品货号: MA02546

第三版

免疫原:Fusion Protein of RAB5A

亚型:IgG2b

克隆:1B6A5

描述:Rab5a belongs to the Rab GTPases family proteins, which play roles in endocytosis, exocytosis, and vesicular transportation. Rab family proteins were recently reported to play important roles in human cancer progression. Rab5a was reported to be overexpressed in breast and ovarian cancers and associated with proliferation and invasion.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF