



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RXRG Monoclonal Antibody

单克隆抗体

产品货号: MA02548

第三版

免疫原:Fusion Protein of RXRG

亚型:IgG2b

克隆:1E7E5

描述:RXRG,also named as,Retinoic acid receptor RXR-gamma,is a 463 amino acid protein,which localize in the nucleus and may form homodimer. Retinoic acid receptors bind as heterodimers to their target response elements in response to their ligands,all-trans or 9-cis retinoic acid,and regulate gene expression in various biological processes. The RAR/RXR heterodimers bind to the retinoic acid response elements (RARE) composed of tandem 5'-AGGTCA-3' sites known as DR1-DR5.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF