



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat SFRP2 Monoclonal Antibody

单克隆抗体

产品货号: MA02551

第三版

免疫原:Fusion Protein of SFRP2

亚型:IgG1

克隆:1F3A1

描述:Secreted frizzled-related protein 2 (Sfrp2) is a secreted glycoprotein molecule containing an N terminal cysteine-rich domain,which is typically 30-50% identical to the putative Wnt-binding site of the frizzled receptor,and a C terminal heparin-binding domain with weak homology to netrins. It has been implicated in diverse cellular processes,including embryogenesis,regulation of cell apoptosis,and cell differentiation. Methylation of this gene is a potential marker for the presence of colorectal cancer.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF