

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLCO1B3 Monoclonal Antibody

单克隆抗体

产品货号: MA02554

第三版

免疫原:Fusion Protein of SLCO1B3

亚型:IgG1

克隆:1D9A4

描述:SLCO1B3,also known as OATP8 or OATP1B3,is a member of Organic Anion Transporting Polypeptides (OATPs) which are sodium-independent organic anion transporters mediating the uptake of a wide range of structurally diverse endogenous and exogenous compounds including bile acids,hormone conjugates,peptides,toxins,as well as a multitude of therapeutic drugs. SLCO1B3 is specifically expressed in liver under normal condition,while its expression has also been observed in cancer tissues like colon,prostate and pancreas. The molecular mass of SLCO1B3 is highly dependent on the glycosylation modification (deglycosylated form of 65 kDa,fully-glycosylated form of 120 kDa).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF