
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SOD2 Monoclonal Antibody

单克隆抗体

产品货号: MA02558

第三版

免疫原:Fusion Protein of SOD2

亚型:IgG2a

克隆:3A6C2

描述:SOD2(superoxide dismutase 2,mitochondrial) is also named as IPOB,MNSOD,SODM,Mn-SOD and belongs to the iron/manganese superoxide dismutase family. It is a marker of mitochondria,which is restricted to the perinuclear area in a cell with aggregate formation of mutant SOD1. It functions as the first line of antioxidant defense against highly reactive superoxide radicals and it appears to be early predictors for survival in septic patients with with MIF. It has 2 isoforms with the molecular weight of 25 kDa and 21 kDa.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF