



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SRD5A1 Monoclonal Antibody

单克隆抗体

产品货号: MA02561

第三版

免疫原:Fusion Protein of SRD5A1

亚型:IgG1

克隆:1D5D8

描述:Steroid 5-alpha-reductase (SRD5A1) catalyzes the conversion of testosterone into the more potent androgen, dihydrotestosterone (DHT). SRD5A1 is a minor component of the reductase activity in prostate although the gene was originally cloned from prostate. On the other hand, SRD5A1 appears to be the predominant isozyme of steroid 5-alpha-reductase in the scalp and elsewhere in the skin.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF