



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human TBL2 Monoclonal Antibody

单克隆抗体

产品货号: MA02567

第三版

免疫原:Fusion Protein of TBL2

亚型:IgA

克隆:1D5C2

描述:TBL2 belongs to a member of the beta-transducin protein family. This gene is commonly deleted in Williams-Beuren syndrome, a developmental disorder with multi-system manifestations caused by haploinsufficiency for contiguous genes deleted in chromosome region 7q11.23. TBL2 encodes a protein with four putative WD40 repeats, which may be involved in some intracellular signaling pathway.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF