

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat TGOLN2 Monoclonal Antibody

单克隆抗体

产品货号: MA02570

第三版

免疫原:Fusion Protein of TGOLN2/TGN46

亚型:IgG1

克隆:1F6D5

描述:TGN46 (TGOLN2),the human homologue of rat Tgn38,is a transmembrane glycoprotein predominantly localized to the TGN (trans-Golgi network). TGN is a major secretory pathway sorting station for proteins and lipids. TGN46 may be involved in regulating membrane traffic to and from TGN. Alternatively spliced transcript variants encode different TGN46 isoforms. TGN46 has an apparent molecular mass of 100-150 kDa,suggesting extensive O- and N-glycosylations.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF